BWISE

EN

INSTRUCTION FOR USE
BWISE Gelling Fiber Dressing

Description
BWISE Gelling Fiber Dressing is a soft, conformable, non-woven pad or ribbon dressing composed of carboxymethyl cellulose fiber and ES-fiber (HDPE/PET). This conformable and highly absorbent dressing absorbs wound fluid and creates a soft gel, maintains a moist environment which supports the body's healing process and aids in the removal of unnecessary material from the wound (autolytic debridement), without damaging newly formed tissue. This primary dressing should be used with a secondary cover dressing.

Indications
BWISE Gelling Fiber Dressing is used for the management of moderately to heavily exuding acute or chronic wounds. It can also be pre-moistened for effective use on dry or lightly draining wounds, including:

• Lower leg ulcers, pressure ulcers (Stage II-IV) and diabetic ulcers;
• Surgical wounds (post-operative, donor sites, wounds left to heal by secondary intent, dermatological);
• Partial thickness burns;
• Traumatic wounds (abrasions and lacerations);
• Oncology wound (if moderately or heavily exuding, superficial or deep).

Contraindication
Patients with allergies to the sodium carboxymethyl cellulose (CMC) should refrain from using this product.

Adverse Reactions
No adverse effects have been reported to date with correct use.

Directions for use
1. Before applying the dressing, cleanse the wound area with an appropriate wound cleaner.

2. Apply the dressing onto the wound and cover with a moisture retentive dressing, like foam dressing, gauze, or other appropriate dressing. BWISE Gelling Fiber Dressing should overlap 1.0 cm or 1/2 inch onto the skin surrounding the wound.
3. When using BWISE ribbon in deep cavity wounds, Only fill deep wounds to 80%, as the ribbon will expand to fill the wound space on contact with wound fluid.
4. In donor sites, BWISE Gelling Fiber Dressing may be left in place for up to 14 days.

5. All wounds should be inspected frequently. Remove the BWISE Gelling Fiber Dressing when clinically indicated (i.e., leakage, excessive bleeding, suspicion of infection).
6. The BWISE Gelling Fiber Dressing is designed to remain in place no t exceed 7 days. The dressing should be changed when it is saturated with wound fluid or if the cover dressing is leaking or the cover dressing's edges are bunching or rolling up.
7. When using BWISE Gelling Fiber Dressing to wounds with lower exudate levels, moisten it with sterile physiological saline solution (0.9%) before changing the dressing. Any gel residue on the wound should be removed when cleansing the wound. Deep wounds in particular should be well irrigated.

Warning
In case of wound swelling or redness for more than 48 hours after use of BWISE Gelling Fiber Dressing, contact a doctor.
Caution
This product is single use only and is supplied sterile. Discard any unused portion of the product after dressing the wound.

If the inner package is damaged, do not use.
Packages and Storage
BWISE Gelling Fiber Dressing is sterile and packed individually. Different sizes are available.
Keep BWISE Gelling Fiber Dressing in a cool dry place for maximum shelf life. The preferred storage temperature should be between 5°C and 35°C.

Sterile
BWISE Gelling Fiber Dressing is sterilized in the original packaging by electron beam radiation.

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